

Innovations...

Vial2Bag™

Vial2Bag needleless reconstitution for IV systems

The Vial2Bag needleless system enables safe and convenient reconstitution and transfer of a drug between a vial or syringe and any solution-filled IV bag. The needle-free system connects to standard IV bags through the IV set port so it can be used with all manufacturers' bags.

The drug vial connects to the Vial2Bag system with a standard vial adapter. The diluent is transferred from the IV bag to the vial for drug reconstitution. The reconstituted drug is then transferred from the vial back to the bag. The system helps provide transfer of the full dose to the IV bag for administration. Standard syringes can also be connected to the Vial2Bag system.

Vial2Bag system attributes

- Fits standard IV systems
- Helps ensure optimum drug utilization
- Safe, needleless reconstitution and transfer
- Secure connection and transfer between a vial or syringe and the IV bag
- Multi-purpose repeat access Luer connection
 - Sterility is maintained on the inside
 - Contact surface is swabable
- Molded from medical-grade materials

Vial2Bag configurations

- The system is delivered sterile in a blister pack.
- The IV spike is covered with a removable protector.
- A tamper-evident feature provides access to connect the IV set spike.
- The system can be kitted with the drug vial and a vial adapter.

Vial2Bag is ideal for

- Hospital pharmacy drug preparation
- Out-patient health care centers
- Clinical uses

World-class technologies from Medimop

The Vial2Bag system was designed by Medimop Medical Projects Ltd., a West Pharmaceutical Services company. Medimop is the world market leader in the design and manufacture of transfer, mixing and administration systems for injectable pharmaceuticals.



Innovations in Reconstitution, Mixing and Transfer Systems

West Pharmaceutical Services' drug reconstitution, mixing and transfer systems enhance the delivery of drugs through advanced technologies.

West's systems are easy to use, help promote safety and compliance and do not require changes to drug manufacturing processes.

To find out how your drug product can benefit from the Vial2Bag system, contact your West Pharmaceutical Services' account representative.

In North America,
call 908-918-1015.

In Europe, call +45 8227-3029.

You may also visit our web site at
westReconstitution.com.

